# 3D Bond+™

Bone Graft Cement

## Instructions for use

(EN) English

(FR) Français

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(IT) Italiano

(DE) Deutsch

(ES) Español

(RO) Română

(PT) Português

(PL) Polski

(HU) Magyar

(DA) Dansk

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### English

#### INSTRUCTION FOR USE

Please read this entire circular before performing the procedure.

This device is for sale by, or on the order of a physician, or a licensed practitioner. The intended users are clinicians who are familiar with augmentation procedures in the dental and maxillofacial field, for example: dental and maxillofacial surgeons

#### COMPOSITION

Each syringe contains 0.5cc medical grade Biphasic Calcium Sulfate and 0.3 mL physiological saline (0.9% sodium chloride for injection).

#### PRODUCT DESCRIPTION

3D Bond+™ is a synthetic osteoconductive, bioresorbable bone graft cement. composed of biphasic Calcium Sulfate in a granulated powder form. 3D Bond+™ function as a scaffold for bone regeneration in dental procedures.

#### INTENDED PURPOSE

3D Bond+™ is intended to fill, augment, or reconstruct bony defects in the oral and maxillofacial region.

#### INDICATION FOR USE

3D Bond+™ is indicated for use in the following ways: by itself in bone regenerative techniques or mixed with other suitable bone filling agents, to prevent particle migration in an osseous defect, and to provide a resorbable barrier over other bone graft material.

#### CLINICAL BENEFITS

3D Bond+™ is osteoconductive and bioactive, allowing a fast-healing time and bone regeneration. Furthermore, the product prevents infiltration of epithelial connective cells and enables soft tissue proliferation above its surface, therefore, no membrane is required, and primary closure is not mandatory.

Use of 3D Bond+™ as graft material results in minimal inflammation while superb handling and high predictability is given.

#### PATIENT PROFILE

3D Bond+™ can be used among males and females mostly from 18-85 years old with diverse background and different socioeconomic status. Smokers and medically compromised patients are treated as well based on risk benefit evaluation for any surgical procedure. The decision and final discretion are made by the clinician.

#### STERILIZATION

3D Bond+™ sterilization is carried out using Gamma irradiation.

#### CONTRAINDICATIONS

- The customary contraindications in oral and maxillofacial surgery with other implant materials should be observed.
- Lack of adequate training of the practitioner is a major risk for the success of the implant procedure.
- Acute and chronic active infection at the site of the implant.
- Serious disturbance of bone metabolism
- Serious bone diseases of endocrine etiology
- Sever or difficult to control diabetes mellitus.
- Immunosuppressive and radiation therapy
- Ongoing treatment with glucocorticoids mineralocorticoid and with agent effecting calcium metabolism
- Malignancies
- Lactating and pregnancy
- Effect on pediatric patients is not known.

#### PRECAUTIONS AND WARNINGS

- The material is designed for single-use only. DO NOT RE-STERILIZE. 3D Bond+ must no longer be used in case of partially opened or defective primary packaging (syringe) or secondary packaging (peel-off-blister) since the sterility of the material is no longer ensured.
- Do not use when the temperature of the product is below 10° c (50° F). If used in low temperatures, wait until the product restores room temperature low temperature will slow the setting reaction of the material).
- It is not recommended to mix 3D Bond+™ with blood.
- 3D Bond+™ has short-term space maintaining characteristics (4-10 weeks). Therefore, it is not indicated to augment large bony defects, or defect when a long-term space maintaining is required.
- The expiration date is printed on the peel-off blister and on the external package. Do not use after indicated expiry date.
- 3D Bond+™ does not possess sufficient mechanical strength to support load bearing defects prior to tissue in-growth. In case load support is required, standard internal and external stabilization technique must be followed to obtain rigid stabilization in all planes.
- 3D Bond+™ must not be used to stabilize screw placement.
- 3D Bond+™ is intended for use by clinicians familiar with bone grafting procedures.
- Possible complications are the same as to be expected of autogenous bone grafting procedures. These may include superficial wound infection, deep wound infection, deep wound infection with osteomyelitis, delayed union, loss of reduction, failure of fusion, loss of bone graft, graft protrusion and /or dislodgment, and general complications that may arise from anesthesia and/or surgery.

- Complication specific to oral/dental use are those as may be typically observed for similar bone grafting procedures and may include, tooth sensitivity, gingival recession, flap sloughing, resorption or ankylosis of the treated root, abscess formation.

### STORAGE

Store at temperature between 5°C (41°F) to 30° C (86°F), avoid contact with a source of heat. Do not store the product in direct sunlight.

### HANDLING

3D Bond+™ is available as a granulated powder packed within a specific double compartment syringe. The syringe contains the biphasic Calcium Sulfate powder in one compartment and a sterile standard saline (0.9% sodium chloride for injection) in the second compartment.

### IMPORTANT

- Read all steps of the instructions before using 3D Bond+™
- It is highly recommended to practice the use of 3D Bond+™ before first use.

### INSTRUCTIONS FOR SITE PREPARATION

- Reflect the mucoperiosteal flap, if required.
- Remove the undesired soft tissue from the exposed bone surface.
- Prepare the defect area for augmentation procedure.

Note that a membrane is not required in most surgical procedures. In addition, primary flap closure is recommended but not essential for proper healing (since soft tissue may grow on top of 3D Bond+™).

### STEP BY STEP

Read all instruction steps (1-3) before using 3D Bond+™

It is recommended that the augmented site will be completely debrided and prepared before activating the material, the material should be injected to the site **immediately** after its activation.

### PREPARATION STEP

Before applying 3D Bond+™ make sure you have the items described below:

- a- A sterile 3D Bond+™ syringe
- b- Dry sterile gauze pads
- c- Sterile periosteal elevator or spatula

### STEP 1.

Place your index finger firmly on syringe cap and slowly push the shaft towards the line marked on the syringe tube until the first plunger reaches the line. This will activate the material and prepare it for ejection.

Note: while pushing the shaft, mild pressure is required.

After completion of above step, **wait 5 seconds** for complete hydration of the cement before removing the syringe cap.

### STEP 2

Remove the cap by twisting and pulling it out.


### STEP 3

Eject the material into the augmented site (socket)

Note: the paste should be in direct contact with the bone and slightly overfilled.

After injecting the material into the site:  
place a dry sterile gauze pad on the material and firmly condense the cement by exerting pressure with a finger above the gauze for 3 seconds.

Replace the finger pressure by applying additional compaction using periosteal elevator or spatula on the gauze for 3-5 second. **The material must be well compacted at the cervical zone of the socket.**

Remove the gauze, remove excess material and slightly shape if required. Press on the material once again with a dry sterile gauze

Remove the gauze and proceed with soft tissue coverage and wound closure.

#### Notes:

- In cases of 4 bony walls socket it is not required to reach the apex zone with the graft cement material, however, the graft cement should be well compacted at the socket cervical level.
- In cases when buccal or lingual plate are missing, it is recommended priorly to reflect a minimal flap to enable the introduction of the graft cement material into the entire cavity of the socket, and compact it firmly in all three dimensional aspects.
- In cases when the graft cement is left exposed to the oral cavity ,in order to prevent the material and volume lose during the healing phase, the exposed surface of the graft must be protected by collagen sponge, or adhesive wound dressing band that should be placed and secured by suturing above the graft cement and the soft tissue for 7-14 days until soft tissue proliferation will cover the exposed area.

#### POST TREATMENT CARE

- Do not place any removeable provisional appliance above the grafted site during the entire healing stage.
- Instructions that include an appropriate analgesic, antibiotic, and home care regimen should be delivered to the patient after surgical intervention.
- When 3D Bond+™ is used for filling bony defects pre implant placement, let the site heal for 3 to 6 months before dental implant placement.
- Prior to implant placement the grafted site must be evaluated to ensure that adequate bone healing has occurred.
- Place the dental implant and abutment according to the cleared indications and instructions for dental implant/abutment system used.

#### RESORPTION PROFILE

3D Bond+™ completely resorbs in approximately 4-10 weeks.

#### WASTE DISPOSAL AFTER USAGE

the disposal of the syringe after use should be in accordance with contaminated medical waste disposal instructions.

## INCIDENT REPORTING

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP shall be made available upon request to [info@augmabio.com](mailto:info@augmabio.com)

The SSCP can be found in the European database on medical devices (Eudamed) at

<https://ec.europa.eu/tools/eudamed>.

Basic UDI-DI of device: 7290014838ABM5143YW

## KEY TO USED CODES

\* ISO 15223-1: Medical Devices - Symbols to be used with medical device labels, labelling and information to be supplied. – Part 1: General requirements

|  |                                                                                                                                                                                                              |
|--|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Manufacturer<br>(5.1.1)*                                                                                                                                                                                     |
|  | <b>Authorized representative in the European Union</b><br>MedNet EC-REP III GmbH, Borkstrasse 10, 48163 Münster, Germany<br>(5.1.2)*                                                                         |
|  | <b>CE Mark</b> with Notified Body identification number.                                                                                                                                                     |
|  | Date and Country of manufacturer<br>(5.1.3, 5.1.11)*                                                                                                                                                         |
|  | Use-by date<br>(5.1.4)*                                                                                                                                                                                      |
|  | Batch code<br>(5.1.5)*                                                                                                                                                                                       |
|  | Catalogue/Part number<br>(5.1.6)*                                                                                                                                                                            |
|  | Importer<br>Indicates the entity importing the medical device into the locale. This symbol shall be accompanied by the name and address of the importing entity adjacent to the symbol.<br>(5.1.8)*          |
|  | Distributor<br>Indicates the entity distributing the medical device into the locale. This symbol shall be accompanied by the name and address of the distributing entity adjacent to the symbol.<br>(5.1.9)* |
|  | Sterilized using irradiation<br>(5.2.4)*                                                                                                                                                                     |
|  | Do not resterilize<br>(5.2.6)*                                                                                                                                                                               |

|                    |                                                                                                                                                                                                                                                    |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                    | Do Not Use if package is damaged<br>(5.2.8)*                                                                                                                                                                                                       |
|                    | Single sterile barrier system with protective packaging inside<br>(5.2.13)*                                                                                                                                                                        |
|                    | Single sterile barrier system with protective packaging outside<br>(5.2.14)*                                                                                                                                                                       |
|                    | Keep away from sunlight<br>(5.3.2)*                                                                                                                                                                                                                |
|                    | Temperature limit<br>Indicates the temperature limits to which the medical device can be safely exposed.<br>(5.3.7)*                                                                                                                               |
|                    | Do not re-use/ single use/ use only once<br>(5.4.2)*                                                                                                                                                                                               |
|                    | Consult electronic instructions for use<br>(5.4.3)*                                                                                                                                                                                                |
|                    | Caution<br>Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that, for a variety of reasons, cannot be presented on the medical device itself.<br>(5.4.4)* |
|                    | Medical Device<br>(5.7.7)*                                                                                                                                                                                                                         |
|                    | Unique Device Identifier<br>(5.7.10)*                                                                                                                                                                                                              |
|                    | Prescription only<br>Caution:<br>US federal law restricts this product for sale by or on the order of a dentist or physician.<br>(21 CFR 801.15 (c) (1)(i) (F), 801.109)                                                                           |
|                    | <b>Authorised representative for Switzerland</b><br>MedNet SWISS GmbH, D4 Platz 4, 6039 Root D4, Switzerland<br>(MedDO)                                                                                                                            |
| Qty                | Quantity                                                                                                                                                                                                                                           |
| 1 Syringe          | 1 Syringe                                                                                                                                                                                                                                          |
| cc                 | Dosage                                                                                                                                                                                                                                             |
| Bone Graft Cement  | Bone Graft Cement                                                                                                                                                                                                                                  |
| 3D Bond+           | 3D Bond+                                                                                                                                                                                                                                           |
| Augma Biomaterials | Augma Biomaterials                                                                                                                                                                                                                                 |

# 3D Bond+™

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0344

**Augma Biomaterials Ltd.**

Alon Hatavor 20 St., P.O.Box 3089

Caesarea Southern Industrial Park

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