3D BOND™

Graft Binder Cement

(EN) English

(FR) Français

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(IT) Italiano

(DE) Deutsch

(ES) Español

(RO) Română

(PL) Polski

(HU) Magyar

(DA) Dansk

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Augma Biomaterials Ltd.

Alon Hatavor 20 St., P.O.Box 3089

Caesarea Southern Industrial Park

3088900, Israel

www.augmabio.com

3D Bond™

Graft Binder Cement

Augma Biomaterials Ltd.
Alon Hatavor 20 St., P.O. Box 3089
Caesarea Southern Industrial Park
3088900, Israel
www.augmabio.com

English

INSTRUCTIONS FOR USE

Please read this entire circular before performing the procedure.

This device is for sale by, or on the order of a physician, or licensed practitioner. The intended users are clinicians who are familiar with augmentation procedures in the dental and maxillofacial field, for example: dental and maxillofacial surgeons

COMPOSITION

Each driver contains 1cc medical grade Biphasic Calcium Sulfate.

PRODUCT DESCRIPTION

3D Bond™ is a synthetic, osteoconductive, bioresorbable graft binder cement composed of biphasic Calcium Sulfate in a granulated powder form. 3D Bond™ functions as a scaffold for bone regeneration in dental procedures.

INTENDED PURPOSE

3D Bond™ is intended to fill, augment, or reconstruct bony defects in the oral and maxillofacial region.

INDICATIONS FOR USE

3D Bond™ is indicated for use in the following ways: by itself in bone regenerative techniques or mixed with other suitable bone filling agents, to prevent particle migration in an osseous defect, and to provide a resorbable barrier over other bone graft material.

CLINICAL BENEFITS

3D Bond™ is osteoconductive and bioactive, allowing a fast-healing time and bone regeneration. Furthermore, the product prevents infiltration of epithelial connective cells and enables soft tissue proliferation above its surface, therefore, no membrane is required, and primary closure is not mandatory. Use of 3D Bond™ as graft material results in minimal inflammation while superb handling and high clinical predictability is given.

PATIENT PROFILE

3D Bond™ can be used among males and females mostly from 18-85 years old with diverse background and different socioeconomic status. Smokers and medically compromised patients are treated as well based on risk benefit evaluation for any surgical procedure. The decision and final discretion are made by the clinician.

STERILIZATION

3D Bond™ sterilization is carried out using γ – irradiation.

CONTRAINDICATIONS

- The customary contraindications in oral and maxillofacial surgery with other implant materials should be observed
- Lack of adequate training of the practitioner is a major risk for the success of the implant procedure.
- Acute and chronic active infection at the site of the implant
- Serious disturbance of bone metabolism
- Serious bone diseases of endocrine etiology
- Severe or difficult to control diabetes mellitus
- Immunosuppressive and radiation therapy
- Ongoing treatment with glucocorticoids mineralocorticoids and with agents affecting calcium metabolism
- Malignancies
- Lactating and pregnancy
- Effect on pediatric patients is not known

PRECAUTIONS AND WARNINGS

- The material is supplied for single use only. DO NOT RE-STERILIZE. 3D Bond™ must no longer be used in case of partially opened or defective primary packaging (driver) or secondary packaging (peel-off-blister) since the sterility of the material is no longer ensured.
- Do not use when the temperature of the product or the saline is below 10°C (50°F). If used in low temperatures, wait until the product restores room temperature (low temperatures will slow the setting reaction of the material).
- It is not recommended to mix 3D Bond™ with blood.
- 3D Bond™ has short-term space maintaining characteristics (4-10 weeks). Therefore, by itself it is not indicated to augment large bony defects, or defects when long-term space maintaining is required.
- In large defects, 3D Bond™ should be mixed with slow reabsorbing granular augmentation materials in a volumetric ratio of 2:1 (for example, mix 1cc 3D Bond™ with 0.5cc other granular material). Other mixing ratios could compromise the graft hardening in-situ. The ratio calculation should be based on the same units (volume or weight).
- The expiration date is printed on the peel-off blister and on the external package. Do not use after indicated expiry date.
- 3D Bond™ does not possess sufficient mechanical strength to support load bearing defects prior to tissue in-growth. In case load support is required, standard internal or external stabilization techniques must be followed to obtain rigid stabilization in all planes.
- 3D Bond™ must not be used to stabilize screw placement.
- 3D Bond™ is intended for use by clinicians familiar with bone grafting procedures.

- Possible complications are the same as to be expected of autogenous bone grafting procedures. These may include: superficial wound infection, deep wound infection, deep wound infection with osteomyelitis, delayed union, loss of reduction, failure of fusion, loss of bone graft, graft protrusion and/or dislodgement, and general complications that may arise from anesthesia and/or surgery.
- Complications specific to oral/dental use are those as may be typically observed for similar bone grafting procedures and may include: tooth sensitivity, gingival recession, flap sloughing, resorption or ankylosis of the treated root, abscess formation.

STORAGE

Store at temperatures between 5°C (41°F) to 30°C (86°F). Avoid contact with a source of heat. Do not store the product in direct sunlight.

HANDLING

3D Bond™ is available as granulated powder packed within a driver. A sterile standard saline (0.9% Sodium Chloride for injection) is required on-site as a setting agent.

IMPORTANT

- Read all steps of the instructions (1-7) before using 3D Bond™.
- It is highly recommended to practice the use of 3D Bond™ before first usage.

INSTRUCTION FOR SITE PREPARATION

- Reflect the mucoperiosteal flap, if required.
- Remove the undesired soft tissue from the exposed bone surface.
- Prepare the defect area for augmentation procedure.

Note that a membrane is not required in most surgical procedures. In addition, primary flap closure is recommended but not essential for proper healing (since soft tissue may grow on top of 3D Bond™). 3D Bond™ can be used also as a membrane above other augmentation materials. In that case, follow the same steps and apply 3D Bond™ over the selected augmentation material in a thickness which is not less than 1mm.

3D Bond™ can be used by mixing with other suitable bone filling agents to prevent particle migration in an osseous defect. Mixing has to be done in a ratio of 2:1 (For example, mix 1cc 3D Bond™ with 0.5cc other granular material). The preferred granule size should not exceed 1.5mm/0.06”.

OPTIONAL METHOD (the purpose of the optional method is to reduce bleeding):

- If there is an accumulation or abundant of bleeding at the surgical site, apply 3D Bond™ in two layers.
- The first layer should be thin. Inject a small amount of 3D Bond™ into the defect area and compress it for few seconds with a sterile, dry gauze pad (a significant decrease of bleeding at the surgical site is expected).
- Apply a second and final 3D Bond™ layer above the first one by compressing it using a dry gauze for few seconds (the second layer is applied in order to fill the defect cavity entirely).
- You may remove excess material as required using a suitable instrument moistened with saline. After removal, press firmly again with the gauze for 1-2 seconds.
- Proceed with soft tissue coverage and wound closure.

STEP BY STEP

Using 3D Bond™ as a Composite Graft

Read all instruction steps (1-7) before using 3D Bond™. It is recommended that the time elapsed from the introduction of saline into 3D Bond™ should not exceed 3 minutes; otherwise, the material will lose its pliability.

PREPARATION STEP

Before applying 3D Bond™, make sure you have the items described:

A - A sterile 3D Bond™ driver

B - Sterile dry gauze pads


C - Sterile medical syringe with sterile saline (Sodium Chloride 0.9% for injection)

D - Sterile medical needle 21-27G (1.25” -1.5” / 31.75mm-38.1mm)

E - Sterile granules of your preferable bone graft (the preferred granules size should not exceed 1.5mm/0.06”)

F - Sterile dish

G - Sterile spatula

STEP 1.

Add 3D Bond™ to the granules within the dish. Use a volume ratio of 2 parts 3D Bond™ to 1 part granules. (For example, add 1cc 3D Bond™ to 0.5cc granules).

STEP 2.

Use a spatula and mix the dry blend thoroughly.

**STEP 3.**

Add 0.8 cc saline for one portion of 1.5cc (granules/ 3D Bond™) for complete wettings of the blend (This is important in order to achieve a complete moistening of the blend). (A)

Use a spatula to mix the blend thoroughly until obtaining a homogeneous creamy texture (The texture is ready when it's well mixed and has a shiny and glossy appearance). (B)


**STEP 4.**

Absorb water surpluses by applying a dry gauze pad to the mixture. Do not use pressure (The texture is ready when it has a matt appearance).

**STEP 5.**


Use the empty 3D Bond™ driver as a carrier for the putty by retracting the piston to line 1 marked on the driver and re-load it with the putty.

Note: Make sure that the time elapsed from the introduction of saline into 3D Bond™ until it is use, will not exceed 3 minutes.

**STEP 6.**

Apply the 3D Bond™ composite putty into the required site by pushing the piston and ejecting the putty. (A)

The composite paste should be in direct contact with the bone and slightly overfilled. (B)

STEP 7.

After injecting the material into the site:

Place a dry, sterile gauze pad on the material and firmly condense the cement by exerting pressure with a finger above the gauze for 3 seconds. Next, replace the finger pressure by applying additional compaction using periosteal elevator or spatula on the gauze for 3-5 second. **The material must be well compacted from all directions.**

Remove the gauze, remove excess material and slightly shape, if required. Press on the material once again with a dry, sterile gauze for 1-2 seconds.

Remove the gauze and proceed with soft tissue coverage and wound closure.

STEP BY STEP

Using 3D Bond™ by Itself

Read all instruction steps (1-7) before using 3D Bond™. It is recommended that the time elapsed from the introduction of saline into 3D Bond™ should not exceed 3 minutes; otherwise, the material will lose its pliability.

3D Bond™ by itself has short-term space maintaining characteristics (4-10 weeks), therefore, by itself is indicated to augment relatively small bone deficiency. For large bony defects, or defects where long-term space maintaining is required, 3D Bond should be used as a composite graft.

PREPARATION STEP

Before applying 3D Bond™, make sure you have the items described:


A - A sterile 3D Bond™ driver

B - Sterile saline (Sodium Chloride 0.9% for injection)

C - Sterile dry gauze pads

D - Sterile medical syringe

E - Sterile medical needle 21G-27G (1.25" -1.5" / 31.75mm-38.1mm)

Please note: Lines 1&2 marked on the driver do not indicate quantities or Measurements. They represent the handling steps.

STEP 1.


Push the piston to line 1 marked on the driver. Twist and pull to remove the driver's cap.

STEP 2.

Insert the syringe's needle through the driver's head and eject the saline from the syringe into the 3D Bond™ driver.

- Eject until excess liquid drops out.
- Do not block the driver's aperture.

**STEP 3.**


Hold the driver's head handles firmly and push the piston towards line 2 to expel excess liquid.

**STEP 4.**

Hold the driver horizontally and remove the driver's head by twisting and pulling.

**STEP 5.**

Dab the driver aperture on dry gauze for 1-2 seconds. This will absorb liquid surplus to obtain the desired viscosity.

**STEP 6.**

Apply the 3D Bond™ paste into the required site by pushing the piston and ejecting the putty. (A)
The paste should be in direct contact with the bone and slightly overfilled. (B)

STEP 7.

After injecting the material into the site:

Place a dry, sterile gauze pad on the material and firmly condense the cement by exerting pressure with a finger above the gauze for 3 seconds. Next replace the finger pressure by applying additional compaction using periosteal elevator or spatula on the gauze for 3-5 second. **The material must be well compacted from all directions.**

Remove the gauze, remove excess material and slightly shape, if required. Press on the material once again for 1-2 seconds with a dry sterile gauze.

Remove the gauze and proceed with soft tissue coverage and wound closure.

POST TREATMENT CARE

- Do not place any removeable provisional appliance above the grafted site during the entire healing stage.
- Instructions that include an appropriate analgesic, antibiotic, and home care regimen should be delivered to the patient after surgical intervention.
- When 3D Bond™ is used for filling bony defects pre implant placement, let the site heal for 3 to 6 months before dental implant placement.
- Prior to implant placement the grafted site must be evaluated to ensure that adequate bone healing has occurred.
- Place the dental implant and abutment according to the cleared indications and instructions for dental implant/abutment system used.

RESORPTION PROFILE

3D Bond™ by itself completely resorbs in approximately 4-10 weeks.

WASTE DISPOSAL AFTER USAGE

The disposal of the syringe after use should be in accordance with contaminated medical waste disposal instructions.

INCIDENT REPORTING

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The SSCP shall be made available upon request to info@augmabio.com

The SSCP can be found in the European database on medical devices (Eudamed) at

<https://ec.europa.eu/tools/eudamed>.

Basic UDI-DI of device: 7290014838ABM5143YW

KEY TO USED CODES

* ISO 15223-1: Medical Devices - Symbols to be used with medical device labels, labelling and information to be supplied. – Part 1: General requirements

	Manufacturer (5.1.1)*
	Authorized representative in the European Union MedNet EC-REP III GmbH, Borkstrasse 10, 48163 Münster, Germany (5.1.2)*
	CE Mark with Notified Body identification number.
	Date and Country of manufacturer (5.1.3, 5.1.11)*
	Use-by date (5.1.4)*
	Batch code (5.1.5)*
	Catalogue/Part number (5.1.6)*
	Importer Indicates the entity importing the medical device into the locale. This symbol shall be accompanied by the name and address of the importing entity adjacent to the symbol. (5.1.8)*
	Distributor Indicates the entity distributing the medical device into the locale. This symbol shall be accompanied by the name and address of the distributing entity adjacent to the symbol. (5.1.9)*
	Sterilized using irradiation (5.2.4)*
	Do not resterilize (5.2.6)*
	Do Not Use if package is damaged (5.2.8)*
	Single sterile barrier system with protective packaging inside (5.2.13)*

	Single sterile barrier system with protective packaging outside (5.2.14)*
	Keep away from sunlight (5.3.2)*
	Temperature limit Indicates the temperature limits to which the medical device can be safely exposed. (5.3.7)*
	Do not re-use/ single use/ use only once (5.4.2)*
	Consult electronic instructions for use (5.4.3)*
	Caution Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that, for a variety of reasons, cannot be presented on the medical device itself. (5.4.4)*
	Medical Device (5.7.7)*
	Unique Device Identifier (5.7.10)*
R only	Prescription only Caution: US federal law restricts this product for sale by or on the order of a dentist or physician. (21 CFR 801.15 (c) (1)(i) (F), 801.109)
Qty	Quantity
2 drivers	2 drivers
cc	Dosage
Graft binder cement	Graft binder cement
3D Bond	3D Bond
Augma Biomaterials	Augma Biomaterials

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